

Your Access To Excellent Results

VAPR[®] S50

Radiofrequency Electrode

The VAPR S50 Electrode offers outstanding access and versatility for soft tissue ablation. The VAPR S50 is a small-diameter suction electrode with an ergonomic integrated handpiece. With its controlled default setting, the VAPR S50 is well suited for meniscal sculpting and sealing, minimizing the jagged appearance that may result from the use of manual instruments.

The VAPR S50 also offers the versatility of precise, fast ablation for ACL debridement. The small diameter and angled tip provide for excellent access to tight spaces, while the robust suction removes bubbles and small debris from the field of vision.

ACCESS & VERSATILITY



ACCESS

- **12% Smaller profile** than the leading competitor's electrode^{1,2}
- **Operative tip angled** for excellent access to tight anatomy

VERSATILITY

- **Controlled default setting** for sealing and sculpting of meniscal tissue
- The power you need for precise, **fast ablation of soft tissue**

1. IMS Data, published June 2006. 2. Measurement Report, D. Archila, March 1, 2007. Data on File.



Excellent Access

Provides excellent access to tighter parts of the joint.

- Low profile tip is 12% smaller than the leading competitor's electrode!¹
- Angled electrode allows smooth access to tight spots.

Versatility & Control

Provides versatile size, angle, and speed for a variety of procedures.

- Controlled default setting for meniscal sculpting, sealing and shaping.
- Independent control of ablation wattage allows for sufficient power for ACL debridement, notchplasty and general tissue ablation.
- 53% faster tissue ablation at max setting than the leading competitor's electrode.³

Exceptional Visualization

VAPR S50 provides robust suction for outstanding visualization.

- Fluid inlet hole designed to manage suction capacity to meet procedural requirements.
- Offers a large-diameter suction tube.

Safety

Advanced electronics provide patient safety and decreased surgical interruption.

- "Seamless Fault Safety Feature" shuts off electrode when it approaches metal and automatically reactivates when safely away.

Support

Provided by representatives with procedural knowledge and the comprehensive product portfolio required during repair.

- Work with trusted, informed professionals.
- Call center support available 24 hours a day, 7 days a week.
- 24-hour console replacement service if and when needed.

ORDERING INFORMATION

Part Number	Description
<i>Featured Electrode</i>	
227355	VAPR S50 Electrode
<i>Other Select VAPR Electrodes</i>	
225370	VAPR S90 Electrode
22730	VAPR Hook Electrode
	w/ Integrated Handpiece
227301	VAPR Side Effect Electrode
	w/ Integrated Handpiece