

*The Performance You Deserve, The Quality You Demand*

# QUICKANCHOR<sup>®</sup> PLUS FAMILY


*NOW fully loaded with ORTHOCORD<sup>™</sup> Suture*



Super QUICKANCHOR Plus Dual Suture



Rotator Cuff QUICKANCHOR Plus Dual Suture



GII QUICKANCHOR Plus

**ORTHOCORD<sup>™</sup>**  
**JUST TIE IT**

- **Strong**  
55 lbs. of tensile strength<sup>1</sup> and 30 lbs. of knot strength<sup>1</sup>
- **45% less stiff<sup>2</sup>**  
In your hands and for your patient
- **Less mass postoperatively**  
Permanent suture with PDS component
- **Experience**  
Developed by ETHICON, INC., with 100 years of suture experience

<sup>1</sup>Data on file

<sup>2</sup>Results of Stiffness Analysis Test TM007-03 performed by ETHICON, INC. CPC. Data on file

 **DePuy Mitek**  
a Johnson & Johnson company

IT'S ABOUT RESULTS



Providing Procedural Solutions...

# QUICKANCHOR<sup>®</sup> PLUS FAMILY

*NOW fully loaded with ORTHOCORD<sup>™</sup> Suture*



## Super QUICKANCHOR Plus Dual Suture

- **Shoulder**  
Rotator cuff repair, capsular repair, biceps repair, deltoid repair
- **Ankle**  
Achilles tendon repair/reconstruction
- **Knee**  
Extra capsular repairs; reattachment of: medial collateral ligament, lateral collateral ligament, posterior oblique ligament, or joint capsule to tibia and joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions

## Rotator Cuff QUICKANCHOR Plus Dual Suture

- **Shoulder**  
Rotator cuff repair

## GII QUICKANCHOR Plus

- **Shoulder**  
Bankart repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/capsulo-labral reconstruction, biceps tenodesis, deltoid repair
- **Ankle**  
Lateral instability, medial instability, Achilles tendon repair/reconstruction, mid-foot reconstruction
- **Foot**  
Hallux valgus reconstruction
- **Wrist**  
Scapholunate ligament
- **Hand**  
Ulnar or lateral collateral ligament reconstruction
- **Elbow**  
Tennis elbow repair, biceps tendon reattachment
- **Knee**  
Extra capsular repairs; reattachment of: medial collateral ligament, lateral collateral ligament, posterior oblique ligament, or joint capsule to tibia and joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions

### ORDERING INFORMATION

Part Number	Description	Order Qty.
222983	GII QUICKANCHOR Plus w/ ORTHOCORD	
222984	Super QUICKANCHOR Plus w/ Dual-Suture ORTHOCORD	
222979	RC QUICKANCHOR Plus w/ Dual-Suture ORTHOCORD	
211093	2.4mm GII Drill Bit	
211502	2.4mm Arthroscopic GII Drill Bit	
211113	2.9mm Super/RC Drill Bit	
211513	2.9mm Arthroscopic Super/RC Drill Bit	
213207	GII/Super/RC Slotted Drill Guide	
214562	Arthroscopic Drill Guide - Saw Tooth	
214646	CordCutter	

Surgeon's Signature

Date

For more information, call your DePuy Mitek representative at 1-800-382-4682 or visit us at [www.depuymitek.com](http://www.depuymitek.com). DePuy Mitek, Inc., 325 Paramount Drive, Raynham, MA 02767

© DePuy Mitek, Inc. 2005. All rights reserved. Printed in the USA. PIN 900857, 11/05 Dyneema<sup>®</sup> Purity is a registered trademark of Royal DSM N.V.

**DePuy Mitek**  
a Johnson & Johnson company

